

## Product Datasheet

### TGF beta Receptor 2 peptide (orb374847)

|                            |                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | TGF beta Receptor 2 peptide                                                                                         |
| <b>Conjugation</b>         | Unconjugated                                                                                                        |
| <b>Tested Applications</b> | HPLC                                                                                                                |
| <b>Target</b>              | TGF beta Receptor 2                                                                                                 |
| <b>Form/Appearance</b>     | Lyophilized                                                                                                         |
| <b>Storage</b>             | Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                                               |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb229827 (TGF beta Receptor 2 antibody)                              |
| <b>Source</b>              | Synthesis                                                                                                           |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                      |