

Product Datasheet

GCN2 (phospho-Thr898) peptide (orb374855)

Description	GCN2 (phospho-Thr898) peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	GCN2 (phospho-Thr898)
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb251471 GCN2 (phospho-Thr898) antibody)
Source	Synthesis
Expiration Date	6 months from date of receipt.