

Product Datasheet

HMGCR (phospho-Ser872) peptide (orb374856)

Description	HMGCR (phospho-Ser872) peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	HMGCR (phospho-Ser872)
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb251472 (HMGCR (phospho-Ser872) antibody)
Source	Synthesis
Expiration Date	6 months from date of receipt.