

Product Datasheet

IL11 peptide (orb374860)

Description	IL11 peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	IL11
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb251476 (IL11 antibody)
Source	Synthesis
Expiration Date	6 months from date of receipt.