




Product Datasheet

AKT antibody (orb37496)

Description	Rabbit polyclonal antibody to AKT
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Synthetic Peptide
Target	AKT1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55686
Uniprot ID	P31749
NCBI	NP_001014432.1 , NP_005154.2 , NP_001014431.1
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MDA-MB435...



Western blot analysis of in MCF-7 cell l...