

Product Datasheet

Caspase 3 peptide (orb374962)

Description	Caspase 3 peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	Caspase 3
Form/Appearance	Lyophilized
Storage	Upon reconstitution, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb340700 (Caspase 3 antibody)
Source	Synthesis
Expiration Date	6 months from date of receipt.