



Product Datasheet

PCTAIRE1 antibody (orb37530)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PCTK1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51091
Uniprot ID	Q9BRL4
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cell lysate...