

Product Datasheet

C1orf112 antibody (FITC) (orb37576)

Description	C1orf112 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Uncharacterized protein C1orf112 protein (301-530AA)
Target	C1orf112
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NSG2
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.