




Product Datasheet

SPHK1 antibody (orb37630)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SPHK1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	42518
Uniprot ID	Q9NYA1
NCBI	NP_001136073.1 , NP_068807.2 , NP_001136074.1 , NP_892010.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of in 293 cell lin...



Western blot analysis of SPHK1 antibody