

Product Datasheet

TNFSF15 Antibody / VEGI / TL1A (orb376640)

Description

Vascular endothelial growth inhibitor (VEGI), also


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P

Applications

Immunohistochemical
staining of FFPE
hum...

Immunogen


A full length human recombinant protein was used as the immunogen for this TL1A antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

Storage

Store the TL1A antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).



Immunohistochemical
staining of FFPE
hum...

Note

For research use only

Application notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the TL1A antibody to be titrated up or down for optimal performance.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

[O95150](#)

Hazard
Information

This TL1A antibody is available for research use only.

Entrez

[9966](#)

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml

Expiration Date

12 months from date of receipt.