



Product Datasheet

Ankyrin R Antibody (orb376692)

Description	Mouse Anti-Human Ankyrin R Monoclonal IgG2B
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Fusion protein amino acids 1-1881 (full-length) of human Ankyrin-R
Target	Ankyrin R
Preservatives	PBS pH 7.4, 50% glycerol, 0.1% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	A 1:100 dilution of SMC-487 was sufficient for detection of Ankyrin R in 20 µg of mouse brain lysate by ECL immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG2B
Clonality	Monoclonal
MW	200kDa
Uniprot ID	P16157
NCBI	NP_000028.3
Entrez	286
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of neuroblas...



Western Blot analysis of rat brain using...