


Product Datasheet

GRP170 Antibody (orb376699)

Description	Mouse Anti-Human GRP170 Monoclonal IgG2bK
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Recombinant Full length GRP170 Protein
Target	GRP170
Preservatives	PBS pH 7.4, 50% glycerol, 0.1% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	1 µg/ml of SMC-232 was sufficient for detection of GRP170 in 20 µg of HEK293 lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG2b
Clonality	Monoclonal
MW	170kDa
Uniprot ID	Q9Y4L1
NCBI	NP_006389.3
Entrez	10525
Dilution Range	WB (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of neuroblas...