



Product Datasheet

HO-1 Antibody (orb376701)

Description	Mouse Anti-Rat HO-1 Monoclonal IgG1 Kappa
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	His-tagged Rat HO-1
Target	HO-1
Preservatives	PBS, 50% glycerol, 0.1% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	1 µg/ml of SMC-234 was sufficient for detection of HO-1 in 10 µg of rat kidney lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1 Kappa
Clonality	Monoclonal
MW	32kDa
Uniprot ID	P06762
NCBI	NP_036712.1
Entrez	24451
Dilution Range	WB (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of fibroblas...



Western Blot analysis of human, mouse, r...