


Product Datasheet


PP5 Antibody (orb376709)

Description	Mouse Anti-Human PP5 Monoclonal IgG1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Full length human PP5 protein
Target	PP5
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	1 µg/ml of SMC-244 was sufficient for detection of PP5 in 15 µg of human A431 lysates by ECL analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	58kDa
Uniprot ID	P53041
NCBI	NP_006238.1
Entrez	5536
Dilution Range	WB (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of cervical

...



Immunofluorescence analysis of HEK293 us...



Western Blot analysis of human A431, HEK...