

Product Datasheet

**beta 4 Spectrin antibody [Out of stock]
(orb376749)**

Description

Mouse monoclonal antibody against

488.

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

ATTO-488

Tested

ICC, IHC, WB

Applications
Immunogen

Fusion protein amino acids 1621-1832 (C-terminal repeats 14 to 15) of human Beta4-spectrin

Target

beta 4 Spectrin

Preservatives

PBS pH 7.4, 50% glycerol, 0.1% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

A 1:100 dilution of SMC-484 was sufficient for detection of Beta 4 Spectrin in 20 µg of mouse brain lysate by ECL immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

MW

200kDa

Uniprot ID

[Q9H254](#)

NCBI

[NP_066022.2](#)

Entrez

[57731](#)

Dilution Range

WB (1:1000)

Expiration Date

12 months from date of receipt.