



Product Datasheet

UVRAG Antibody (orb377178)


Description	Rabbit Anti-Human UVRAG Polyclonal
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Synthetic peptide from the mid-protein of human UVRAG
Target	UVRAG
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-606 was sufficient for detection of UVRAG on HeLa cell lysates using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	78.1kDa
Uniprot ID	Q9P2Y5
NCBI	NP_003360.2
Entrez	7405
Dilution Range	WB (1:1000); ICC/IF (1:100)
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of neuroblas...



Western Blot analysis of human hela_293t...



Western Blot analysis of rat liver using...