



Product Datasheet

GABARAP Antibody (orb377268)


Description	Rabbit Anti-Human GABARAP Polyclonal
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Synthetic peptide from the N-terminal of human GABARAP
Target	GABARAP
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-620 was sufficient for detection of GABARAP on mouse kidney lysates using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	14-16kDa
Uniprot ID	O95166
NCBI	NP_009209.1
Entrez	11337
Dilution Range	WB (1:1000); ICC/IF (1:100)
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of neuroblas...



Western Blot analysis of human 293t usin...



Western Blot analysis of mouse kidney us...