

Product Datasheet

p53 antibody [Out of stock] (orb377398)

Description

Rabbit polyclonal antibody against p53

Species/Host	Rabbit
Reactivity	Human
Conjugation	ATTO-594
Tested Applications	ICC, WB
Immunogen	Synthetic peptide from the C-terminal of human tumor protein p53.
Target	p53
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-682 was sufficient for detection of p53 on HeLa cell lysates using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	53kDa
Uniprot ID	P04637
NCBI	NP_000537.3
Entrez	7157
Dilution Range	WB (1:1000); ICC/IF (1:100)
Expiration Date	12 months from date of receipt.