

Product Datasheet

BDNF antibody [Out of stock] (orb377467)

Description

Rabbit polyclonal antibody against BDNF

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	ATTO-390
Tested Applications	ICC, WB
Immunogen	Synthetic peptide from the N-terminal of human BDNF
Target	BDNF
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-703 was sufficient for detection of BDNF on 293T Rapamycin-treated lysates using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	27.9kDa
Uniprot ID	P23560
NCBI	NP_001137277.1
Entrez	627
Dilution Range	WB (1:1000); ICC/IF (1:100)
Expiration Date	12 months from date of receipt.