

Product Datasheet

Choline Acetyltransferase antibody (orb377500)

Description

Rabbit polyclonal antibody against

to RPE.

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

RPE

Tested

ICC, WB

Applications
Immunogen

Synthetic peptide from the N-terminal to the mid-protein of human Choline O-Acetyltransferase

Target

Choline Acetyltransferase

Preservatives

PBS, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

A 1:1000 dilution of SPC-706 was sufficient for detection of Choline Acetyltransferase on mouse brain lysates using Goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

82.5kDa

Uniprot ID

[P28329](#)

NCBI

[NP_001136401.1](#)

Entrez

[1103](#)

Dilution Range

WB (1:1000); ICC/IF (1:100)

Expiration Date

12 months from date of receipt.