

Product Datasheet

PIK-293 (orb377542)

Description	PIK-293 is a PI3K inhibitor, mostly for PI3Kδ with IC50 of 0.24 μM, 500-, 100- and 50-fold less po...
Note	For research use only
Formula	C ₂₂ H ₁₉ N ₇ O
MW	397.43
CAS Number	[900185-01-5]
SMILES	CC1=C2C(=CC=C1)N=C(N(C2=O)C3=CC=CC=C3C)CN4C5=C(C=N4)C(=NC=N5)N
Expiration Date	12 months from date of receipt.