

Product Datasheet

PLG antibody (orb37767)

Description	Rabbit polyclonal antibody to PLG
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PLG
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	90569
Uniprot ID	P00747
NCBI	NP_000292.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

