

Product Datasheet

Gp2 antibody (orb37774)

Description	Gp2 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Mouse Pancreatic secretory granule membrane major glycoprotein GP2 protein (22-531AA)
Target	Gp2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9D733
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.