

Product Datasheet

Fibrinogen gamma antibody (orb37784)

Description Rabbit polyclonal antibody to Fibrinogen gamma

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IF, WB

Immunogen Synthetic Peptide

Target FGG

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 51512

Uniprot ID [P02679](#)


NCBI [NP_000500.2](#), [NP_068656.2](#)

Dilution Range WB: 1:1000, FACS: 1:10-50, IF/ICC: 1:10-50

Expiration Date 12 months from date of receipt.



Immunofluorescence analysis of HepG2 cel...



Western blot analysis of A431 cell line ...