



Product Datasheet

FERMT3 antibody (orb377840)

Description	Rabbit polyclonal antibody to FERMT3.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	Fermt3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	75635
Uniprot ID	Q8K1B8
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse spleen ly...