



Product Datasheet

PACSIN2 antibody (orb377854)

Description	Mouse monoclonal antibody to PACSIN2.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	PACSIN2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	55739
Uniprot ID	Q9UNF0
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane1) Daudi w...