

Product Datasheet

PPP1R3E antibody (orb377859)

Description	Rabbit polyclonal antibody to PPP1R3E.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	PPP1R3E
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	30644
Uniprot ID	Q9H7J1
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

