


Product Datasheet

CA12 antibody (Biotin) (orb377870)

Description	Goat polyclonal antibody against CA12 conjugated to Biotin.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	carbonic anhydrase XII (aa188-199), Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Shipped at ambient temperature, store at -20°C
Note	For research use only
Application notes	WB: Approx 45kDa band observed in Human kidney lysates (calculated MW of 39.5kDa according to NP_001209.1). Recommended concentration: 0.1-0.3µg/ml.
Protein Sequence	CSHLQHVYKKGQE.
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human Kidney ly...