


Product Datasheet

PTGR2 antibody (orb377873)

Description	Goat polyclonal antibody to PTGR2.
Species/Host	Goat
Reactivity	Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PTGR2 / ZADH1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 38kDa band observed in Pig Testis lysates (calculated MW of 38.4kDa according to Pig XP_001929109.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	PDNINEGQVQVRT
Clonality	Polyclonal
MW	38.5
Entrez	145482
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Pig Testis lysa...