

Product Datasheet

NLRP7 antibody (orb377874)

Description	Goat polyclonal antibody to NLRP7.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	NLRP7
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments in Human Ovary and Uterus lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml).
Protein Sequence	QGDIDNFHDDVT
Clonality	Polyclonal
MW	115; 112; 118
Entrez	199713
Dilution Range	ELISA: 1:32000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.