


Product Datasheet

CTDSP1 antibody (Biotin) (orb377876)

Description	Goat polyclonal antibody against CTDSP1 conjugated to Biotin.
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CTDSP1, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Shipped at ambient temperature, store at -20°C
Note	For research use only
Application notes	WB: Approx 30kDa band observed in Mouse Testis lysates (calculated MW of 29.3kDa according to Mouse NP_694728.1). Recommended concentration: 0.1-0.3µg/ml.
Protein Sequence	EAKAQSDKIC.
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Testis ly...