


Product Datasheet

TPO antibody (Biotin) (orb377879)

Description	Goat polyclonal antibody against TPO conjugated to Biotin.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	thyroid peroxidase, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Shipped at ambient temperature, store at -20°C
Note	For research use only
Application notes	WB: Approx 110kDa band observed in Human Thyroid Gland lysates (calculated MW of 103kDa according to NP_000538.3). Recommended concentration: 1-3µg/ml.
Protein Sequence	ELRCGKHQAVGTSPQ.
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human Thyroid G...