


Product Datasheet

CD274 antibody (Biotin) (orb377881)

Description	Goat polyclonal antibody against CD274 conjugated to Biotin.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CD274 / PD-L1, Biotinylated
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Shipped at ambient temperature, store at -20°C
Note	For research use only
Application notes	WB: Approx 40kDa band observed in Human Heart lysates (calculated MW of 33.3kDa according to NP_054862.1). Recommended concentration: 0.3-1µg/ml.
Protein Sequence	CKKQSDTHLEET.
Dilution Range	ELISA: 1:128000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human Heart lys...