




Product Datasheet

PCSK9 antibody (orb37797)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	PCSK9
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	74286
Uniprot ID	Q8NBP7
NCBI	NP_777596.2
Dilution Range	WB: 1:1000, IHC-P: 1:10-50, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of PCSK9 antibody ...



Flow cytometric analysis of PCSK9 antibo...