

Product Datasheet

SLC25A18 antibody (HRP) (orb37801)

Description	SLC25A18 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Mitochondrial glutamate carrier 2 protein (246-287AA)
Target	SLC25A18
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9H1K4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.