


Product Datasheet

DDX58 antibody (orb378028)

Description	Rabbit polyclonal antibody to DDX58.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide of human DDX58
Target	DDX58
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	095786
Entrez	23586
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:20-100
Expiration Date	12 months from date of receipt.



Western blot analysis of DDX58 expression...