



Product Datasheet

HAVCR2 antibody (orb378097)


Description	Rabbit polyclonal antibody to HAVCR2.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	KLH-conjugated synthetic peptide of human TIM3
Target	HAVCR2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q8TDQ0, Q8VIM0
Entrez	171285, 84868
Dilution Range	WB: 1:500-1000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of TIM3 expression...



Immunohistochemical analysis of TIM3 sta...



Immunofluorescent analysis of TIM3 stain...