



Product Datasheet

HNRNPR antibody (orb378106)

Description	Rabbit polyclonal antibody to HNRNPR.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide of human hnRNP R
Target	HNRNPR
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	O43390
Entrez	10236
Dilution Range	WB: 1:500-1000, IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.




Western blot analysis of hnRNP R express...



Immunohistochemical analysis of hnRNP R

...



Immunofluorescent analysis of hnRNP R st...