




Product Datasheet

SMYD3 antibody (orb37812)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	SMYD3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	49097
Uniprot ID	Q9H7B4
NCBI	NP_001161212.1 , NP_073580.1
Dilution Range	WB 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of in CEM cell lin...



Western blot analysis of HT-29 cell line...