



Product Datasheet

MAP1LC3A antibody (orb378164)


Description	Rabbit polyclonal antibody to MAP1LC3A.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	Recombinant full length protein of human LC3A/B
Target	MAP1LC3A; MAP1LC3B
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q6XVN8 , Q91VR7 , Q9GZQ8 , Q9H492 , Q9CQV6 , Q62625
Entrez	64862 , 66734 , 81631 , 67443 , 84557 , 362245
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF/ICC: 1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of LC3A/B expressi...



Immunohistochemical analysis of LC3A/B s...



Immunofluorescent analysis of LC3A/B sta...