



Product Datasheet

MME antibody (orb378180)

Description	Rabbit polyclonal antibody to MME.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	Recombinant full length protein of human CD10
Target	MME
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q61391 , P07861 , P08473
Entrez	24590 , 17380 , 4311
Dilution Range	WB: 1:500-2000, IP: 1:20-50, FACS: 1:20-50
Expiration Date	12 months from date of receipt.



Western blot analysis of CD10 expression...



Immunohistochemical analysis of CD10 sta...