



Product Datasheet

CD3E antibody (orb378203)


Description	CD3E antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CD3E. Please contact us for the exact immunogen sequence. The peptide is available as orb12965.
Target	CD3E
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	23 kDa
Dilution Range	IHC-P: 1:50-300
Expiration Date	12 months from date of receipt.



IHC-P staining of rat lung tissue using ...



IHC-P image of rat lung tissue using CD3...



WB analysis of rat lymph node (lane 1), ...