



Product Datasheet

Tubulin beta 3 antibody (orb378207)


Description	Tubulin beta 3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Tubulin beta 3. Please contact us for the exact immunogen sequence. The peptide is available as orb12969.
Target	Tubulin beta 3
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	50 kDa
Dilution Range	IF/ICC: 1:100-800, IHC-P: 1:100-800, WB: 1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse brain (la...



IHC-P image of rat spinal cord tissue us...



WB analysis of mouse brain (lane 1), rat...