


Product Datasheet

Casein Kinase Iggamma2 antibody (orb378251)

Description	Rabbit polyclonal antibody to Casein Kinase Iggamma2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:1000 IHC-P: 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...