


## Product Datasheet

### MRCK beta antibody (orb378332)

|                            |                                         |
|----------------------------|-----------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to MRCK beta |
| <b>Species/Host</b>        | Rabbit                                  |
| <b>Reactivity</b>          | Human                                   |
| <b>Conjugation</b>         | Unconjugated                            |
| <b>Tested Applications</b> | IHC-P, WB                               |
| <b>Note</b>                | For research use only                   |
| <b>Clonality</b>           | Polyclonal                              |
| <b>Dilution Range</b>      | WB: 1:1000 IHC-P: 1:200                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.         |



Western blot analysis of COLO205 whole c...