


## Product Datasheet

### USP36 antibody (orb378439)

|                            |                                     |
|----------------------------|-------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to USP36 |
| <b>Species/Host</b>        | Rabbit                              |
| <b>Reactivity</b>          | Human, Mouse                        |
| <b>Conjugation</b>         | Unconjugated                        |
| <b>Tested Applications</b> | IHC-P, WB                           |
| <b>Note</b>                | For research use only               |
| <b>Clonality</b>           | Polyclonal                          |
| <b>Dilution Range</b>      | WB: 1:1000 IHC-P: 1:200             |
| <b>Expiration Date</b>     | 12 months from date of receipt.     |



Western blot analysis of HuvEc whole cel...