


## Product Datasheet

### ZNF262 antibody (orb378451)

|                            |                                      |
|----------------------------|--------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to ZNF262 |
| <b>Species/Host</b>        | Rabbit                               |
| <b>Reactivity</b>          | Human, Mouse                         |
| <b>Conjugation</b>         | Unconjugated                         |
| <b>Tested Applications</b> | IHC-P, WB                            |
| <b>Note</b>                | For research use only                |
| <b>Clonality</b>           | Polyclonal                           |
| <b>Dilution Range</b>      | WB: 1:1000 IHC-P: 1:200              |
| <b>Expiration Date</b>     | 12 months from date of receipt.      |



Western blot analysis of HepG2 whole cel...