



Product Datasheet

GIT1 antibody (orb37915)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	GIT1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	84341
Uniprot ID	Q9Y2X7
NCBI	NP_054749.2 , NP_001078923.1
Dilution Range	WB 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse heart tissue.