



Product Datasheet

DNMT3L Antibody (orb379608)


Description	Mouse monoclonal antibody to DNMT3L
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant fragment of human DNMT3L (AA: 147-386) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	43.6kDa
Entrez	29947
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Line graph illustrates the Ag-Ab r...



Western blot analysis of HEK293 cell lys...



Western blot analysis of recombinant pro...