



Product Datasheet

DDX58 Antibody (orb379676)


Description	Mouse monoclonal antibody to DDX58
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC
Immunogen	Purified recombinant fragment of human DDX58 (AA: 789-925) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgM
Clonality	Monoclonal
MW	106.6kDa
Entrez	23586
Dilution Range	ICC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Line graph illustrates about the Ag-Ab r...



Flow cytometric analysis of HeLa cells u...



Immunofluorescence analysis of MCF-7 cel...